

### PRE-OWNED OCT'S

Refurbished, like new, patient ready!

### What Other Doctors Say

*Luis has helped me on several occasions in both repairing and purchasing equipment. He has been accommodating at all times. After purchasing an OCT from Used Eye Equipment, he arranged to have an excellent tech come in and set up with multiple instructional training sessions for my staff. His response time to attend to any questions or issues has been prompt. Very easy to work with and quality equipment and service.*  
~ Dr. Jacob N NY



**Cirrus 5000**

Current Zeiss model, capable of Anterior Segment Premiere and Angiography. Operates exclusively on Windows 7 or 10.



### **Cirrus 500**

Current Zeiss starter OCT, capable of Anterior Segment Premiere. Lacks an LSLO, meaning there is no real-time view of the fundus during alignment. Not capable of Angiography. Operates exclusively on Windows 7.



## **Cirrus 4000**

Original full-featured Cirrus model, with anterior segment capability. It can be found with either Windows 7 or 10 operating software.



## **Cirrus 400**

Original Zeiss starter OCT, capable of anterior segment imaging. Lacks an LSLO, meaning there is no real-time view of the fundus during alignment. It can be found with either Windows XP or 7 operating software.



## **Optovue RTVue-100**

\*Lombart \$18,495 Our Price \$10,950

Original Optovue spectral-domain OCT. Available with anterior segment imaging, and Optovue's Total Corneal Power (TCP) software; not available with i Wellness software. It can be found with either Windows XP or 7 operating software.



## Optovue iVue

Current Optovue spectral-domain OCT model. Capable of anterior segment imaging and Optovue i Wellness. It can be found with either Windows XP or 7 operating software.

Our Price \$15,950



## Optovue iFusion

Current Optovue model combining the iVue OCT and iCam NM fundus camera. Not possible to capture both modalities simultaneously. It can be found with either Windows XP or 7 operating software. Rarely available pre-owned.



## Optovue iScan

Compact and user-friendly, iScan's streamlined design integrates the operator interface, display, and scan head into one compact console. The iScan is available in two configurations: Essential and Comprehensive (comprehensive adding ant. seg. Technology).





## Topcon Triton

Current Topcon OCT/camera combination unit, and the only retinal Swept Source OCT instrument on the US market. Capable of OCT, NM, FA, and FAF photography. Operates on web-based software available on all Windows Platforms including Windows 7 and above.



## Topcon Maestro

Current Topcon model, capable of both NM fundus photography and OCT with anterior segment imaging. Automated and very user-friendly. Rarely available pre-owned.

**For a Bigger Return on Your Investment**

**Call, text, WhatsApp: 949-362-1150**

**[info@usedeyeequipment.com](mailto:info@usedeyeequipment.com)**