



## Topcon KR 8000 Auto Ref/Keratometer

Reconditioned • Calibrated • Patient Ready



The KR 8000 utilizes breakthrough electro-optical and rotary prism technology, the KR-8000

SUPRA is capable of obtaining perfectly accurate keratometry and/or refraction measurements through pupils as narrow as 2.0 mm!

This means easier and more precise diagnostic results when dealing with glaucoma, asymmetric pupils, and elderly patients. The new SUPRA-Series instruments have also been redesigned for even greater ease of operation and patient comfort and are as much as 50% faster than previous generation units. Improved printing speed, auto measurement, and data output capabilities make this the perfect choice in today's fast-paced practice.

**For a Bigger Return on Your Investment**

**Call, text, WhatsApp: 949-362-1150**

**[info@usedeyeequipment.com](mailto:info@usedeyeequipment.com)**