



## **The Oculus Easyfield**

**Refurbished • Calibrated • Patient Ready**

The Oculus Easyfield is a portable, compact automated perimeter for visual field screening and threshold testing: 30-2, 24-2, 10-2, individual points or sectors. Has built-in CCD camera, integrated printer and software for patient data management & storage for up to 40,000 exams in hand control unit. The unit only weighs 15 pounds.

### **Special Features**

- Small, compact and lightweight
- Classic perimeter principles
- Industry standard printout
- Laptop or hand unit control
- 10 to 30 degree test range
- Screening and threshold strategies
- Grayscale, 3D and change maps

**For a Bigger Return on Your Investment**

**Call, text, WhatsApp: 949-362-1150**

**[info@usedeyeequipment.com](mailto:info@usedeyeequipment.com)**

